



**Signatera™**  
Residual disease test (MRD)

Detect residual disease early.  
Treat with confidence.

## MRD Program

Use Signatera after surgery to evaluate the need for adjuvant chemotherapy and avoid unnecessary treatment




## Surveillance Program

Use Signatera alongside CEA to detect recurrence earlier while it may still be resectable, and reduce false positives

\*Medicare draft coverage for MRD Program: Stage II-III colon cancer and Stage IIA rectal cancer patients

\*Medicare draft coverage for Surveillance Program: Stage II-III colorectal cancer patients




- Each patient gets a personalized custom-built assay to reflect their unique tumor
- Signatera targets 16 tumor-derived variants for each patient
- Clonal mutations are selected that occurred early and persist throughout tumor evolution


# MRD Program

Evaluate the need for adjuvant chemotherapy after surgery to avoid unnecessary treatment

Treatment decision based on risk assessment<sup>5</sup>




Signatera status is the only significant predictor of risk




Signatera validated with high PPV and NPV<sup>1</sup>

PPV: Positive Predictive Value  
NPV: Negative Predictive Value



Series of assays in the adjuvant window

MRD Program (first 6 months)



\*Medicare draft coverage for MRD Program: Stage II-III colon cancer and Stage IIA rectal cancer patients

# Surveillance Program

Detect recurrence while it may be resectable, and triage indeterminate nodules


## Significant lead time over CT and CEA

Determine response to treatment earlier than previously possible

**Signatera detects recurrence earlier than CT imaging<sup>1</sup>**

Average lead time vs CT | 8.7 months

Maximum lead time vs CT | 16.5 months



**Signatera is more accurate than CEA with fewer false positives<sup>1,9,10</sup>**

\*Patient-level specificity 98%. Test-level specificity 99.7%.

## Actionable results<sup>5</sup>

**ctDNA ⊕ High risk**

Consider directed imaging (PET/MRI) to locate the disease while potentially resectable

>97% of patients will relapse without treatment

**ctDNA ⊖ Reduced risk**

Continue monitoring with reassurance


12-14% patients will relapse. Patients who remain negative 2 years post treatment have risk reduced to 3%.



Medicare draft coverage\*

Serial testing with same frequency as CEA


**Surveillance Program (after 6 months)**



\*Medicare draft coverage for Surveillance Program: Stage II-III colorectal cancer patients


# Personalized monitoring and MRD assessment across solid tumors

## Signatera validated across tumor types



## How Signatera works

### The personalized and tumor informed approach



# Experts in cell-free DNA. Optimized for oncology.



Reproductive health



Oncology



Organ transplantation

>2M  
cfDNA tests  
performed

>100  
clinicians, PhD's,  
and scientists

>50  
peer-reviewed  
publications

>80  
patients issued  
or allowed

CAP  
accredited

CLIA  
certified

Learn more at [NateraOncology.com](http://NateraOncology.com)

## REFERENCES

1. Reinert T, Henriksen TV, Christensen E, et al. Analysis of plasma cell-free DNA by ultradeep sequencing in patients with stages I to III colorectal cancer. *JAMA Oncol.* 2019;5(8):1124–1131.
2. Coombes RC, Page K, Salari R, et al. Personalized detection of circulating tumor DNA antedates breast cancer metastatic recurrence. *Clin Cancer Res.* 2019;25(14):4255–4266.
3. Abbosh C, Birkbak NJ, Wilson GA, et al. Phylogenetic ctDNA analysis depicts early-stage lung cancer evolution. *Nature.* 2017;545:446–451.
4. Christensen E, Birkenkamp-Demtroder K, Sethi H, et al. Early detection of metastatic relapse and monitoring of therapeutic efficacy by ultra-deep sequencing of plasma cell-free DNA in patients with urothelial bladder carcinoma. *J Clin Oncol.* 2019;37(18):1547–1557.
5. National Comprehensive Cancer Network. Colon Cancer V. Plymouth Meeting, PA: National Comprehensive Cancer Network. 2019.
6. Aoyama, Oba K, Honda M, et al. Impact of postoperative complications on the colorectal cancer survival and recurrence: analyses of pooled individual patients' data from three large phase III randomized trials. *Cancer Med.* 2017;6(7):1573–1580.
7. Yothers G, O'Connell MJ, Lopatin M, et al. Validation of the 12-gene colon cancer recurrence score in NSABP C-07 as a predictor of recurrence in patients with stage II and III colon cancer treated with fluorouracil and leucovorin (FU/LV) and FU/LV plus oxaliplatin. *J Clin Oncol.* 2013;31(36):4512–4519.
8. Sinicrope FA, Foster NR, Thibodeau SN, et al. DNA mismatch repair status and colon cancer recurrence and survival in clinical trials of 5-Fluorouracil-based adjuvant therapy. *J Natl Cancer Inst.* 2011;103(11):863–875.
9. Purandare NC, Dua SG, Arora A, et al. Colorectal cancer — patterns of locoregional recurrence and distant metastases as demonstrated by FDG PET/CT. *Indian J Radiol Imaging.* 2010;20(4):284–288.
10. Advance Staff. Colorectal cancer monitoring. *Elite Healthcare.* 2016;25(9):14. <https://www.elitecme.com/resource-center/laboratory/colorectal-cancer-monitoring/>



201 Industrial Road, Suite 410 | San Carlos, CA 94070 | [nateraoncology.com](http://nateraoncology.com) | 1.650.489.9050

This test was developed by Natera, Inc., a laboratory certified under the Clinical Laboratory Improvement Amendments (CLIA). This test has not been cleared or approved by the US Food and Drug Administration (FDA). Although FDA does not currently clear or approve laboratory-developed tests in the US, certification of the laboratory is required under CLIA to ensure the quality and validity of the tests. © 2019 Natera, Inc. 2010\_12\_06\_NAT-8020037