


# Patients with fever?


Consider TTM therapy,  
using the Arctic Sun™ Temperature  
Management System.

Treatment with the Arctic Sun™ Temperature Management System has been shown to<sup>1</sup>:

- **Reduce percent time febrile.** Patients **spent less time** febrile at only 8% vs 42% with blankets.
- **Reduce fever** burden by 75% vs traditional water-circulating blankets.
- **Reduce time** to attain normothermia. Patients spent more time normothermic at 59% vs only 3% with blankets.

<sup>1</sup>Mayer, SA. Crit Care Med 2004 Vol. 32, No. 12.





**Patient trend indicator** reflects the rate of change in the patient's temperature over the previous five minutes.

**Precise core targeting** monitors patient temperature every second against the target making micro-adjustments 30x per hour.

**Watch it work**  
Active response technology controls and monitors changes of the three core treatment components.

**Therapy control**

- Protocol
- Time elapsed
- Start treatment

#### Indications for use

The Arctic Sun™ Temperature Management System is a thermal regulating system, indicated for monitoring and controlling patient temperature in adult and pediatric patients of all ages.

#### Contraindications

- There are no known contraindications for the use of a non-invasive thermoregulatory system.
- Do not place ArcticGel™ Pads on skin that has signs of ulcerations, burns, hives or rash.
- Do not remove the fabric release liner of the Neonatal ArcticGel™ Pad and expose the hydrogel.
- Do not place ArcticGel™ Pads on immature (non-keratinized) skin or premature babies.
- While there are no known allergies to hydrogel materials, caution should be exercised with any patient with a history of skin allergies or sensitivities.

#### Warnings

- When using the Arctic Sun™ Temperature Management System, note that all other thermal conductive systems, in use while warming or cooling with this device may interfere with patient temperature control.
- The Arctic Sun™ Temperature Management System is not intended for use in the operating room environment.

#### Cautions

- Due to underlying medical or physiological conditions, some patients are more susceptible to skin damage from pressure and heat or cold. Patients at risk include those with poor tissue perfusion or poor skin integrity due to edema, diabetes, peripheral vascular disease, poor nutritional status, steroid use or high dose vasopressor therapy. Examine the patient's skin under the ArcticGel™ Pads.
- Skin injury may occur as a cumulative result of pressure, time and temperature.
- Carefully remove ArcticGel™ Pads from the patient's skin at the completion of use. Aggressive removal or removal of cold pads from the patient's skin may result in skin tears.
- The rate of temperature change and potentially the final achievable patient temperature is affected by many factors. Treatment application, monitoring and results are the responsibility of the attending physician. If the patient does not reach target temperature in a reasonable time or the patient is not able to be maintained at the target temperature, the skin may be exposed to low or high water temperatures for an extended period of time, which may increase the risk for skin injury.

Please consult Operator's Manual for detailed safety information and instructions for use.  
Federal Law (USA) restricts this device to sale by or on the order of a physician.

Medivance, Inc.  
321 South Taylor Ave., Suite 200  
Louisville, Colorado 80027 U.S.  
[arcticsun@bd.com](mailto:arcticsun@bd.com)



**bd.com**

