



Penumbra System®

ELEVATING PERFORMANCE WITH

Penumbra

JET[®] 7

Reperfusion Catheter with Standard Tip

Penumbra

+ ENGINE[®]



Penumbra 

Penumbra
JET 7
Standard Tip



.072" LUMEN

20 TRANSITIONS
for trackability and navigation

ARTICULATING MARKER BAND
designed to improve tip softness

SUPERIOR FLEXIBILITY
enabled by progressive distal coil wind

ENHANCED PUSHABILITY
featuring Quad-Wire technology

FULL LENGTH PTFE LINER
*designed for durability
with adjunctive devices*

2.16 mm




132 cm length

HIGHEST TRF

WITH PENUMBRA JET 7 STANDARD TIP
POWERED BY PENUMBRA ENGINE

$$\text{TRF} = \text{Catheter Tip Area} \times \text{Vacuum Level}$$

Thrombus removal force (TRF) equates to catheter tip force and drives clot capture.



ORDERING INFORMATION

Penumbra System®

Catalog Number	Description	Proximal OD (F) (in.)	Distal OD (mm)	Proximal ID (in.)	Distal ID (in.)	Working Length (cm)
Aspiration Kits						
5MAXJET7BKIT	Penumbra JET® 7 Reperfusion Catheter with Standard Tip + Penumbra Hi-Flow Tubing	6 (.085)	2.16	.072	.072	132
5MAXJETDKIT	Penumbra JET D Reperfusion Catheter + Penumbra Hi-Flow Tubing	6 (.080)	1.65	.064	.054	138
5MAXACE068KIT	ACE™ 68 Reperfusion Catheter + Penumbra Hi-Flow Tubing	6 (.080)	2.03	.068	.068	132
5MAXACE132KIT	ACE60 Reperfusion Catheter + Penumbra Hi-Flow Tubing	6 (.080)	1.80	.068	.060	132
4MAXCKIT	4MAX™ Reperfusion Catheter + Penumbra Hi-Flow Tubing	6 (.080)	1.42	.064	.041	139
3MAXCKIT	3MAX Reperfusion Catheter + Penumbra Hi-Flow Tubing	4.7 (.062)	1.27	.043	.035	160

Reperfusion Catheters		PENUMBRA SYSTEM – Indication for Use Penumbra Reperfusion Catheters and Separators. As part of the PENUMBRA SYSTEM, the Reperfusion Catheters and Separators are indicated for use in the revascularization of patients with acute ischemic stroke secondary to intracranial large vessel occlusive disease (within the internal carotid, middle cerebral – M1 and M2 segments, basilar, and vertebral arteries) within 8 hours of symptom onset. Patients who are ineligible for intravenous tissue plasminogen activator (IV t-PA) or who fail IV t-PA therapy are candidates for treatment. Penumbra 3D REVASCLARIZATION DEVICE As part of the PENUMBRA SYSTEM, the Penumbra 3D REVASCLARIZATION DEVICE is indicated for use in the revascularization of patients with acute ischemic stroke secondary to intracranial large vessel occlusive disease (within the internal carotid, middle cerebral – M1 and M2 segments) within 8 hours of symptom onset. Patients who are ineligible for intravenous tissue plasminogen activator (IV t-PA) or who fail IV t-PA therapy are candidates for treatment. Penumbra Aspiration Tubing As part of the PENUMBRA SYSTEM, the Penumbra Sterile Aspiration Tubing is indicated to connect the Penumbra Reperfusion Catheters to the Penumbra Aspiration Pump. Penumbra Aspiration Pump The Penumbra Aspiration Pump is indicated as a vacuum source for Penumbra Aspiration Systems. Contraindications There are no known contraindications. Warnings The device is intended for single use only. Do not resterilize or reuse. Resterilization and/or Reuse may result in ineffective catheter coating lubrication, which may result in high friction and the inability to access the target neuro vasculature location. Do not use kinked or damaged devices. Do not use open or damaged packages. Return all damaged devices and packaging to the manufacturer/distributor. Do not use automated high-pressure contrast injection equipment with the Penumbra Reperfusion Catheter because it may damage the device. Confirm vessel diameter, and select an appropriate size Penumbra Reperfusion Catheter. Do not use in arteries with diameters smaller or equal to the distal outer diameter of the Penumbra Reperfusion Catheters. Refer to the Reperfusion Catheter labeling for dimensional information. Do not advance, retract or use any component of the PENUMBRA SYSTEM against resistance without careful assessment of the cause using fluoroscopy. If the cause cannot be determined, withdraw the device or system as a unit. Unrestrained torquing or forced insertion of the catheter, revascularization device, or SEPARATOR against resistance may result in damage to the device or vessel. Do not use the PENUMBRA SYSTEM with a pump other than the Penumbra Aspiration Pump. The Penumbra 3D REVASCLARIZATION DEVICE has not been evaluated in patients with angiographic evidence of pre-existing arterial injury. Potential Adverse Events Possible complications include, but are not limited to, the following: acute occlusion; hematoma or hemorrhage at access site; death; intracranial hemorrhage; hemorrhage; infection (at access site); distal embolization; ischemia (cardiac and/or cerebral); embolus (air, foreign body, thrombus, plaque); aneurysm perforation; false aneurysm formation; neurological deficits including stroke; vessel spasm, thrombosis, dissection, perforation or rupture; air embolism; emboli.
5MAXJET7B	Penumbra JET 7 Reperfusion Catheter with Standard Tip	
5MAXJETD	Penumbra JET D Reperfusion Catheter	
5MAXACE068	ACE68 Reperfusion Catheter	
5MAXACE132	ACE60 Reperfusion Catheter	
4MAXC	4MAX Reperfusion Catheter	
3MAXC	3MAX Reperfusion Catheter	
Revascularization Device		
PSR3D	3D Revascularization Device™	
Delivery Microcatheter		
VEL160STR	Velocity® Microcatheter	
Separator™ Devices		
PSF054	5MAX Separator	
PSF041	4MAX Separator	
3MAXS	3MAX Separator	
Aspiration Accessories		
PMXENGN	Penumbra ENGINE®	
PAPS3	Penumbra ENGINE Canister	

BENCHMARK™ BMX™ 96 Access System and Delivery Catheters

Catalog Number	Description	Working Length (cm)
BMX9680BER105	BENCHMARK BMX96 Access System, 80 cm Straight, 105 cm BER 6 F Select™ Catheter	80
BMX9690BER125	BENCHMARK BMX96 Access System, 90 cm Straight, 125 cm BER 6 F Select Catheter	90
BMX9690SIM125	BENCHMARK BMX96 Access System, 90 cm Straight, 125 cm SIM 6 F Select Catheter	90
BMX9690SIMV125	BENCHMARK BMX96 Access System, 90 cm Straight, 125 cm SIM-V 6 F Select Catheter	90
BMX9690MBER125	BENCHMARK BMX96 Access System, 90 cm MP, 125 cm BER 6 F Select Catheter	90
BMX9690MSIM125	BENCHMARK BMX96 Access System, 90 cm MP, 125 cm SIM 6 F Select Catheter	90
BMX9690MSIMV125	BENCHMARK BMX96 Access System, 90 cm MP, 125 cm SIM-V 6 F Select Catheter	90
BMX96100BER125	BENCHMARK BMX96 Access System, 100 cm Straight, 125 cm BER 6 F Select Catheter	100
BMX96100SIM125	BENCHMARK BMX96 Access System, 100 cm Straight, 125 cm SIM 6 F Select Catheter	100
BMX9680	BENCHMARK BMX96 Delivery Catheter 80 cm	80
BMX9690	BENCHMARK BMX96 Delivery Catheter 90 cm	90
BMX96100	BENCHMARK BMX96 Delivery Catheter 100 cm	100






Neuron MAX® 088 6 F Lumen Long Sheaths

Catalog Number	Description	Working Length (cm)
(Crosscut Valve, RHV, and Dilator Included)		
PNML6F088804	Neuron MAX 088 6 F Long Sheath, 80/4 Straight	80
PNML6F088804M	Neuron MAX 088 6 F Long Sheath, 80/4 MP	80
PNML6F088904	Neuron MAX 088 6 F Long Sheath, 90/4 Straight	90
PNML6F088904M	Neuron MAX 088 6 F Long Sheath, 90/4 MP	90
PNML6F0881004	Neuron MAX 088 6 F Long Sheath, 100/4 Straight	100
PNML6F0881004M	Neuron MAX 088 6 F Long Sheath, 100/4 MP	100

6 F Select Catheters

Catalog Number	Description	Working Length (cm)
PNS6F105BER	6 F Select Catheter, 105 BER	105
PNS6F125SIM	6 F Select Catheter, 125 SIM	125
PNS6F125SIMV	6 F Select Catheter, 125 SIM-V	125
PNS6F125BER	6 F Select Catheter, 125 BER	125

Tip Shapes

BENCHMARK BMX96, Neuron MAX 088		Select Catheter
		
Straight	MP	H1
		
		SIM
		
		SIM-V

Caution: Federal (USA) law restricts these devices to sale by or on the order of a physician. Prior to use, please refer to the Instructions for Use for complete product indications, contraindications, warnings, precautions, potential adverse events, and detailed instructions for use. Please contact your local Penumbra representative for more information.

Copyright ©2019–2021 Penumbra, Inc. All rights reserved. The Penumbra P logos, Penumbra System, Penumbra JET, Penumbra ENGINE, ACE, MAX, 3D, 3D Revascularization Device, Separator, Velocity, BENCHMARK, BMX, BMX96, Neuron, Neuron MAX, and Select are registered trademarks or trademarks of Penumbra, Inc. in the USA and other countries. All other trademarks are the property of their respective owners. 13741, Rev. D 04/21 USA

PENUMBRA SYSTEM – Indication for Use
Penumbra Reperfusion Catheters and Separators As part of the PENUMBRA SYSTEM, the Reperfusion Catheters and Separators are indicated for use in the revascularization of patients with acute ischemic stroke secondary to intracranial large vessel occlusive disease (within the internal carotid, middle cerebral – M1 and M2 segments, basilar, and vertebral arteries) within 8 hours of symptom onset. Patients who are ineligible for intravenous tissue plasminogen activator (IV t-PA) or who fail IV t-PA therapy are candidates for treatment. **Penumbra 3D REVASCLARIZATION DEVICE** As part of the PENUMBRA SYSTEM, the Penumbra 3D REVASCLARIZATION DEVICE is indicated for use in the revascularization of patients with acute ischemic stroke secondary to intracranial large vessel occlusive disease (within the internal carotid, middle cerebral – M1 and M2 segments) within 8 hours of symptom onset. Patients who are ineligible for intravenous tissue plasminogen activator (IV t-PA) or who fail IV t-PA therapy are candidates for treatment. **Penumbra Aspiration Tubing** As part of the PENUMBRA SYSTEM, the Penumbra Sterile Aspiration Tubing is indicated to connect the Penumbra Reperfusion Catheters to the Penumbra Aspiration Pump. **Penumbra Aspiration Pump** The Penumbra Aspiration Pump is indicated as a vacuum source for Penumbra Aspiration Systems. **Contraindications** There are no known contraindications. **Warnings** The device is intended for single use only. Do not resterilize or reuse. Resterilization and/or Reuse may result in ineffective catheter coating lubrication, which may result in high friction and the inability to access the target neuro vasculature location. Do not use kinked or damaged devices. Do not use open or damaged packages. Return all damaged devices and packaging to the manufacturer/distributor. Do not use automated high-pressure contrast injection equipment with the Penumbra Reperfusion Catheter because it may damage the device. Confirm vessel diameter, and select an appropriate size Penumbra Reperfusion Catheter. Do not use in arteries with diameters smaller or equal to the distal outer diameter of the Penumbra Reperfusion Catheters. Refer to the Reperfusion Catheter labeling for dimensional information. Do not advance, retract or use any component of the PENUMBRA SYSTEM against resistance without careful assessment of the cause using fluoroscopy. If the cause cannot be determined, withdraw the device or system as a unit. Unrestrained torquing or forced insertion of the catheter, revascularization device, or SEPARATOR against resistance may result in damage to the device or vessel. Do not use the PENUMBRA SYSTEM with a pump other than the Penumbra Aspiration Pump. The Penumbra 3D REVASCLARIZATION DEVICE has not been evaluated in patients with angiographic evidence of pre-existing arterial injury. **Potential Adverse Events** Possible complications include, but are not limited to, the following: acute occlusion; hematoma or hemorrhage at access site; death; intracranial hemorrhage; infection (at access site); distal embolization; ischemia (cardiac and/or cerebral); embolus (air, foreign body, thrombus, plaque); aneurysm perforation; false aneurysm formation; neurological deficits including stroke; vessel spasm, thrombosis, dissection, perforation or rupture; air embolism; emboli.

PENUMBRA ENGINE – Indication for Use
The PENUMBRA ENGINE is indicated as a vacuum source for Penumbra Aspiration Systems. **Contraindications** There are no contraindications. **Warnings** The device is intended for single use only. Do not resterilize or reuse. Resterilization and/or Reuse may result in ineffective catheter coating lubrication, which may result in high friction and the inability to access the target neuro vasculature location and/or may compromise the structural integrity of the device. Do not use kinked or damaged devices. Do not use open or damaged packages. Return all damaged devices and packaging to the manufacturer/distributor. Use prior to the "Use By" date. Use the BENCHMARK BMX96 System in conjunction with fluoroscopic visualization. Do not advance or withdraw the BENCHMARK BMX96 System against resistance without careful assessment of the cause using fluoroscopy. If the cause cannot be determined, withdraw the device. Unrestrained moving or torquing the device against resistance may result in damage to the vessel or device. Maintain a constant infusion of an appropriate flush solution. If flow through the device becomes restricted, do not attempt to clear the lumen by infusion. Remove and replace the device. **Potential Adverse Events** Possible complications include, but are not limited to, the following: acute occlusion; air embolism; death; distal embolization; emboli; false aneurysm formation; hematoma or hemorrhage at puncture site; infection; intracranial hemorrhage; ischemia; neurological deficits including stroke; vessel spasm, thrombosis, dissection, or perforation.

NEURON MAX System – Indication for Use
The NEURON MAX System is indicated for the introduction of interventional devices into the peripheral, coronary, and neuro vasculature. **Contraindications** There are no known contraindications. **Warnings** The NEURON MAX System should only be used by physicians who have received appropriate training in interventional techniques. **Potential Adverse Events** The device is intended for single use only. Do not resterilize or reuse. Resterilization and/or Reuse may result in ineffective catheter coating lubrication, which may result in high friction and the inability to access the target neuro vasculature location and/or may compromise the structural integrity of the device. Do not use kinked or damaged devices. Do not use open or damaged packages. Return all damaged devices and packaging to the manufacturer/distributor. Use prior to the "Use By" date. Use the NEURON MAX System in conjunction with fluoroscopic visualization. Do not advance or withdraw the NEURON MAX System against resistance without careful assessment of the cause using fluoroscopy. If the cause cannot be determined, withdraw the device. Unrestrained moving or torquing the device against resistance may result in damage to the vessel or device. Maintain a constant infusion of an appropriate flush solution. If flow through the device becomes restricted, do not attempt to clear the lumen by infusion. Remove and replace the device. **Potential Adverse Events** Possible complications include, but are not limited to, the following: acute occlusion; air embolism; death; distal embolization; emboli; false aneurysm formation; hematoma or hemorrhage at puncture site; infection; intracranial hemorrhage; ischemia; neurological deficits including stroke; vessel spasm, thrombosis, dissection, or perforation.

Penumbra, Inc. USA
One Penumbra Place
Alameda, CA 94502
USA
1.888.272.4606
T 1.510.748.3200
F 1.510.748.3232
order@penumbrainc.com
info@penumbrainc.com
www.penumbrainc.com

Penumbra 