


HOPE IN LIVER TRANSPLANTATION

COCHRANE REVIEW 2023¹

Superior clinical outcomes compared to cold storage with HOPE¹

HOPE was associated with **improvement in clinically relevant outcomes** in a **meta-analysis of 4 RCTs** evaluating the efficacy and safety of this machine perfusion technique compared with SCS.

HOPE vs. SCS¹



DBD & DCD:
IMPROVED graft survival

HR 0.45; 95% CI 0.23–0.87
P = 0.02



ECD-DBD: **LOWER** number
of serious adverse events

OR 0.45; 95% CI 0.22–0.91
P = 0.03



DCD: **REDUCED** damage
to the bile ducts

OR 0.31; 95% CI 0.11–0.92
P = 0.03

HOPE Key Facts

Simple and effective hypothermic liver perfusion

Perfusate	Belzer MPS® (actively oxygenated)
Timing & Duration	Safe perfusion up to 11.5 hours ² and a total preservation window up to 20 hours ³
Temperature	4–12 °C
Perfusate oxygenation	≥ 60 kPa > 450 mmHg
Technique	◆ Single portal perfusion ◆ Allows for various HOPE protocols

Timepoint:
end-ischemic
in recipient
center



Liver perfusion with the VitaSmart™ System

HOPE CLINICAL RESULTS

Evidence for HOPE in liver transplantation

	N	EAD	Graft Survival	Patient Survival	Complications/SAEs	Biliary complications	NAS
DBD							
Cziganý – ECD-DBD (2021): improved 1-year graft survival + fewer 90-day major complications⁵							
		n.s.	1-year P=0.029	1 year n.s.	90d CD ≥3a P=0.036	All n.s.	
HOPE	23	17%	91%	91%	44%	17%	X
SCS	23	35%	78%	83%	74%	26%	X
Ravaioli – ECD-DBD (2022): lower rate of EAD + improved 1-year graft survival⁶							
		P=0.007	1-year P=0.03		90d CD ≥3b n.s.	All n.s.	
HOPE	55	13%	98%	X	22%	16%	X
SCS	55	35%	87%	X	33%	13%	X
Schlegel (2023): less liver-related graft loss⁷							
			Liver-rel GL P=0.004	Patient-rel GL n.s.	90d CD ≥3a n.s.	AS	
HOPE	85	17%	0%	5%	52%	17%	1%
SCS	85	46%	7%	1%	54%	21%	4%
DCD							
Van Rijn (2021): fewer non-anastomotic biliary strictures⁸							
					90d CD ≥3a	AS	P=0.03
HOPE	78	26%	X	X	N=101	X	6%
SCS	78	40%	X	X	N=132	X	18%
DBD & DCD							
Eden (2024): excellent 5-year patient survival and death-censored graft survival in higher risk donors (ECD-DBD, DCD) following HOPE, similar to benchmark results for unperfused, standard criteria DBDs⁹							
			5 year Death Censored	5 year			24 months
DBD	768	X	90%	85%	X	X	2.50%
DCD	434	X	81%	78%	X	X	11.50%
Reich – (under review): lower rate of EAD + fewer complications¹⁰							
		P<0.001	n.s.	n.s.	1-year CD ≥3a per pt P=0.048	n.s.	Clinically Relevant (DCD)
HOPE	109	20%	95%	97%	2.2	21%	11%
SCS	110	37%	93%	96%	1.9	28%	19%

Additional relevant clinical outcomes



Reduction in acute cellular rejection with HOPE compared to SCS¹¹

All grafts
(N=211)
OR 0.54
95% CI 0.29–1
P=0.05

DCD grafts
(N=153)
OR 0.37
95% CI 0.14–1
P=0.05

Maspero – meta-analysis: 6 studies (2023)¹¹



Less HCC recurrence with HOPE in DCD grafts, compared to unperfused DBD grafts¹²

HCC recurrence
HOPE DCD
6%
vs. SCS DBD 26%; P=0.002

5Y tumor-free survival
HOPE DCD
92%
vs. SCS DBD 73%; P=0.027

Mueller – matched cohort study: N=70 per group (2020)¹²

PROTECTIVE MECHANISMS OF HOPE

Ischemia: mitochondrial injury in liver cells^{13,14}


Mitochondria suffer **severe metabolic changes during ischemia**, which become evident and aggravated after reperfusion (reoxygenation). As a consequence of normothermic reperfusion (on a device or at implantation), the full picture of the **inflammatory ischemia reperfusion injury (IRI)** cascade develops, including release of **reactive oxygen species (ROS)** from mitochondria.¹³

Ischemia (hypoxia)

– O₂

Anaerobic metabolism

- Rapid ATP loss
- Accumulation of succinate + NADH



Reperfusion (reintroduction of O₂)

+ O₂

Re-establishing aerobic metabolism¹⁵

- Slow ATP recovery
- ROS and cytokine release
- Ubiquitous inflammation of the entire organ

Graphic modified from Schlegel A, et al. *J Hepatol.* 2022.¹⁴

HOPE acts before the injury occurs^{14,15,17}

SUPPORTS

NADH + ATP reloading



DECREASES

ROS release & IRI-related inflammation



REDUCES

downstream complications after LT

Viability testing during HOPE with FMN


In addition to ROS, **flavinmononucleotide (FMN)** is a complex-I molecule that is released into the perfusate upon reintroduction of oxygen. The amount of FMN released **increases with the degree of mitochondrial injury** and can be easily measured by spectroscopy, making it a **useful target for viability testing**.^{15,18}



- FMN from complex-I can be **quantified** during HOPE^{15,19}
- Perfusate levels of FMN correspond with clinical outcomes¹⁹
- First assessment after **30 min** of perfusion^{15,19}
- Internationally validated approach published²⁰
- FMN validated as a predictor of EAD and patient survival²¹


BRIDGE TO HOPE US TRIAL RESULTS

Bridge to HOPE randomized trial demonstrates promising benefits for HOPE over SCS alone¹⁰



EAD

HOPE = 20%, vs.
SCS = 37%, $p < 0.001$



Steroid resistant rejection

HOPE = 1%,
SCS = 8%, $p = 0.03$



Overall CD ≥3a complications/ patient

HOPE = 138 complications
in 73 patients,
SCS = 170 complications
in 78 patients, $p = 0.048$



Liver-related CD ≥3a complications/ patient

HOPE = 51 complications
in 36 patients,
SCS = 78 complications
in 54 patients, $p = 0.014$



Median hospital length of stay

HOPE = 8 days,
SCS = 11 days, $p = 0.04$



For more information about HOPE
and other innovative products,
visit www.bridgetolife.com



Download a digital copy of this and
other HOPE and product documents.

Abbreviations

AS: anastomotic strictures; **ATP:** adenosine triphosphate; **CD:** Clavien-Dindo classification; **CI:** confidence interval; **DBD:** donation after brain death; **DCD:** donation after circulatory death; **EAD:** early allograft dysfunction; **ECD:** extended-criteria donor; **FMN:** flavin mononucleotide; **GL:** graft loss; **HCC:** hepatocellular carcinoma; **HOPE:** hypothermic oxygenated machine perfusion; **HR:** Hazard Ratio; **IRI:** ischemia/reperfusion injury; **LT:** liver transplantation; **MPS:** machine perfusion solution; **NADH:** nicotinamide adenine dinucleotide; **NAS:** non-anastomotic strictures; **OR:** odds ratio; **RCT:** randomized controlled trial; **ROS:** reactive oxygen species; **SCS:** static cold storage

References

1. Tingle SJ, et al. Machine perfusion in liver transplantation. *Cochrane Database Syst Rev*. 2023;9(9):CD014685.
2. Boteon A, et al. Eleven hours of hypothermic oxygenated machine perfusion (HOPE) for complex liver retransplantation: A case report. *Artificial Organs*. 2023;00:1-3.
3. Brüggewirth IMA, et al. Prolonged hypothermic machine perfusion enables daytime liver transplantation – an IDEAL stage 2 prospective clinical trial. *EClinicalMedicine*. 2024;68:102411.
4. Czigan Z, et al. Hypothermic Oxygenated Machine Perfusion Reduces Early Allograft Injury and Improves Post-transplant Outcomes in Extended Criteria Donation Liver Transplantation From Donation After Brain Death: Results From a Multicenter Randomized Controlled Trial (HOPE ECD-DBD). *Ann Surg*. 2021;274(5):705-712.
5. Czigan Z, et al. Improved outcomes after hypothermic oxygenated machine perfusion in liver transplantation-Long-term follow-up of a multicenter randomized controlled trial. *Hepatol Commun*. 2024;8(2):e0376.
6. Ravaioli M, et al. Hypothermic oxygenated perfusion in extended criteria donor liver transplantation – A randomized clinical trial. *Am J Transplant*. 2022;22(10):2401-2408.
7. Schlegel A, et al. A multicenter randomized-controlled trial of hypothermic oxygenated perfusion (HOPE) for human liver grafts before transplantation. *J Hepatol*. 2023;78(4):783-793.
8. van Rijn R, et al. Hypothermic Machine Perfusion in Liver Transplantation – A Randomized Trial. *N Engl J Med*. 2021;384(15):1391-1401.
9. Eden J, et al. Long-term outcomes after hypothermic oxygenated perfusion and transplantation of 1,202 donor livers in a real-world setting (HOPE-REAL study). *JHEP*. 2024.06.035.
10. Reich DJ, Mao S, Satish S, Nassar A, et al. Hypothermic Oxygenated Perfusion by End-Ischemic Portal-Venous Approach for Liver Transplantation: Results from the Bridge to HOPE. Multicenter Randomized Controlled Trial (clinicaltrials.gov: NCT05045794). Manuscript under review for publication.
11. Maspero M, et al. Acute rejection after liver transplantation with machine perfusion versus static cold storage: A systematic review and meta-analysis. *Hepatology*. 2023;78(3):835-846.
12. Mueller M, et al. Hypothermic Oxygenated Liver Perfusion (HOPE) Prevents Tumor Recurrence in Liver Transplantation From Donation After Circulatory Death. *Ann Surg*. 2020;272(5):759-765.
13. Chouchani ET, et al. Ischaemic accumulation of succinate controls reperfusion injury through mitochondrial ROS. *Nature*. 2014;515(7527):431-435.
14. Schlegel A, Porte R, Dutkowski P. Protective mechanisms and current clinical evidence of hypothermic oxygenated machine perfusion (HOPE) in preventing post-transplant cholangiopathy. *J Hepatol*. 2022;76(6):1330-1347.
15. Schlegel A, et al. Hypothermic oxygenated perfusion protects from mitochondrial injury before liver transplantation. *EBioMedicine*. 2020;60:103014.
16. Wang L, et al. Flavin Mononucleotide as a Biomarker of Organ Quality-A Pilot Study. *Transplant Direct*. 2020;6(9):e600.
17. Boteon Y, et al. Preventing Tumour Recurrence after Liver Transplantation: The Role of Machine Perfusion. *Int J Mol Sci*. 2020;21(16):5791.
18. Panconesi R, et al. Viability Assessment in Liver Transplantation-What Is the Impact of Dynamic Organ Preservation? *Biomedicines*. 2021;9(2):161.
19. Muller X, et al. Novel Real-time Prediction of Liver Graft Function During Hypothermic Oxygenated Machine Perfusion Before Liver Transplantation. *Ann Surg*. 2019;270(5):783-790.
20. Sun K, et al. Quantifying Flavin mononucleotide: an internationally validated methodological approach for enhanced decision making in organ transplantation. *eBioMedicine*. 2025;116:105761.
21. Dingfelder J, et al. Validation of mitochondrial FMN as a predictor for early allograft dysfunction and patient survival measured during hypothermic oxygenated perfusion. *Liver Transpl*. 2025;31(4):476-488.

VitaSmart™ is CE Marked and available for sale in several markets outside of the United States.

VitaSmart™ has not been cleared by the FDA, and the safety and effectiveness of VitaSmart™ has not been established in the US.