




C-Armor[®] DRAPE

PROTECT PATIENTS. SAFEGUARD STERILITY.

The patented, easy-to-use solution, specifically designed to protect sterility during procedures with multiple C-arm swings.



Help keep the sterile field sterile

Protecting patients from the risk of contamination during surgical procedures is a top priority for today's caregivers, who seek solutions that preserve and maintain sterility while easing workflow.

AORN 2017 recommendations: "Unsterile equipment ... should be covered on the top, bottom, and sides with sterile barrier materials before being introduced to or brought over a sterile field."¹

Improving the sterility of draping techniques can help reduce the risk of increases in hospital costs and post-operative infections.²

Using C-Armor

Effective sterility maintenance

- Adhere the C-Armor drape above the sterile field line to help maintain the sterile field.
- Help reduce the risk of contamination with five-sided coverage of equipment introduced to the sterile field.

Ease of use


- Attach and detach the C-Armor drape simply and quickly during procedures with the fastener tabs.
- Allow for multiple rotations with a single drape due to an innovative design.
- Reposition the drape easily, as the C-arm rotates, with a patented inverted pouch design.

Efficiency

- Help promote procedural standardization by using only one barrier per procedure.
- Reduce the need for excess materials, steps and workflow changes.



C-Armor: Form and Function



1. Strong adhesive ensures a secure grip on the table above the sterile field line.
2. Three fastener tabs offer easy deployment and security when the drape is not deployed.
3. Drape lies flat against the table when not in use for minimal disruption of surgical access.
4. Aluminum band assists in supporting the drape above the sterile field level when deployed.
5. Drape helps protect the C-arm from fluid damage.
6. Five-sided protection helps maintain the sterile field.

Order Information

Code	Description	Quantity
5523	C-Armor	12
5536	C-Armor Kit (Includes C-Armor Cover 5434: 42" x 74")	12

Discover how the C-Armor drape can help you protect patients by preserving sterility. Contact your sales representative or read more at <http://tidiproducts.com>. For U.S. and foreign patent information, visit go.tidiproducts.com/patents

References

1. Guidelines for Perioperative Practice, 2018, AORN, Inc. I.Tech.IV.a.1.
2. Gershkovich, G., Hampton, D., C-Arm Fluoroscopy Draping: A Comparison of Three Techniques, Philadelphia PA



TIDI Products, LLC • 570 Enterprise Drive • Neenah, WI 54956 USA

Customer Service 800.521.1314 • International +1.920.751.4036

Fax 800.837.7770 • www.tidiproducts.com