




THE TRUE HYBRID SOLUTION THAT PROMOTES HEALING

GORE® SYNECOR Biomaterial features Gore core technology to deliver high strength and lasting reliability.

- High-risk ventral repairs: High BMIs, comorbidities and recurrent hernias.
- Versatile, easy-to-use hybrid tri-layer solution.
- Tissue-building Gore 3D PGA:TMC web scaffold promotes rapid vascularity within seven days¹ and tissue ingrowth within one month.*
- Intraperitoneal and preperitoneal configurations.
- Laparoscopic and open ventral hernia repair including robotic procedures.
- For bridging, where there is a need for permanent strength.

MORE FOR YOUR PATIENTS

More size options now available



GORE® SYNECOR Preperitoneal Biomaterial

* Data on file 2015; W.L. Gore & Associates, Inc; Flagstaff, AZ.



Configurations



GORE® SYNECOR Intraperitoneal Biomaterial

Catalogue number	Size
GKFC12	12 cm diameter*
GKFW1015	10 cm × 15 cm†
GKFW1520	15 cm × 20 cm†
GKFR2025	20 cm × 25 cm
GKFR2030	20 cm × 30 cm

Trocar compatibility

GKFR2030 is designed to fit through a 12 mm trocar incision.

GORE® SYNECOR Preperitoneal Biomaterial

Catalogue number	MORE SIZE OPTIONS	Minimum recommended trocar or trocar incision size
GKWC09	9 cm diameter* Now available	8 mm trocar
GKWW1015	10 cm x 15 cm†	10 mm trocar
GKWR1215	12 cm x 15 cm	11 mm trocar
GKWW1520	15 cm x 20 cm†	12 mm trocar
GKWR2025	20 cm x 25 cm	15 mm trocar
GKWR2030	20 cm x 30 cm	15 mm trocar‡
GKWR3030	30 cm x 30 cm Now available	18 mm trocar incision
GKWR3040	30 cm x 40 cm Now available	18 mm trocar incision
GKWR3535	35 cm x 35 cm Now available	21 mm trocar incision
GKWR4040	40 cm x 40 cm Now available	24 mm trocar incision
GKWR4050	40 cm x 50 cm Now available	24 mm trocar incision

* Circle.

† Oval.

‡ Wetting recommended.

1. Crawford N. Assessment of Vascularity via Micro CT in Various Patch Devices. Flagstaff, AZ: W. L. Gore & Associates, Inc.; 2016. [Final study report]. 2344TL.

 Consult Instructions for Use
eifu.goremedical.com

Refer to *Instructions for Use* at eifu.goremedical.com for a complete description of all applicable indications, warnings, precautions and contraindications for the markets where this product is available. FX only

Products listed may not be available in all markets.

GORE, *Together, improving life*, SYNECOR and designs are trademarks of W. L. Gore & Associates.

© 2021, 2022 W. L. Gore & Associates, Inc. 22435535-EN MARCH 2022

W. L. Gore & Associates, Inc.

Flagstaff, AZ 86004

Asia Pacific +65 67332882 **Australia/New Zealand** 1800 680 424 **Europe** 00800 6334 4673

United States 800 437 8181 928 779 2771 goremedical.com

