



Small Bowel Video Capsule Endoscopy: Training and Competency Assessment

System Requirements: PillCam Desktop Reader & Software v9.0

Specification	Minimum	Recommended
CPU	Intel® Core™2 Quad • AMD equivalent	
OS (Operating System)	PillCam Desktop: <ul style="list-style-type: none"> • Windows® 7 SP1 and above, 32 & 64 bit • Windows® 8, 32- & 64-bit • Windows® 8.1, 32- & 64-bit • Windows® 10 	
RAM	1 GB	<ul style="list-style-type: none"> • 3 GB for 32-bit Oss • 6 GB for 64-bit OSs
Disk Space	<ul style="list-style-type: none"> • 1 GB for software • 10 GB available for data 	
Screen (Color and Aspect Ratio / Resolution)	24-bit color: <ul style="list-style-type: none"> • 4:3 / 1024 x 768 • 16:9 / 1366 x 768 • 16:10 (8:5) / 1280 x 800 	32-bit color
Removable Media	DVD-ROM	DVD-RW
USB	USB 2 <ul style="list-style-type: none"> • USB 2 hub (all DR2s and DR3s must connect through USB 2 hub unless Windows® 7 • USB disk or external drivers (for data transfers) must be partitioned as NTFS 	
Software	PDF Reader software such as Adobe® Acrobat	

[IT Guide: Information On Installation & Technical Requirements](#)