

Venue Family



Transducer Guide

The Venue™, Venue Go™ and Venue Fit™ systems provide technologies to help deliver a simple, fast and precise solution to the point of care ultrasound community. Powered by cSound™ Architecture these durable transducers provide exceptional image quality with feature designs that help enhance ease of use, ergonomics and patient comfort.

Transducer	Description	Applications	Footprint	Bandwidth	Biopsy Guide	System Compatibility
 3Sc-RS	Wide-band phased array	Abdomen, Adult Cardiac, Adult Cephalic, Gynecology, Obstetrics, Ophthalmic, Pediatric/Neonatal, Pediatric Cardiac, Thoracic/Pleural, Urology	18.4 x 23.7 mm FOV 120°	Probe Frequency Range: 1.1 – 4.7 MHz	N/A	Venue, Venue Go and Venue Fit
 6S-RS	Wide band phased array	Abdomen, Adult Cardiac, Gynecology, Pediatric/Neonatal, Pediatric Cardiac, Neonatal Cephalic, Thoracic/Pleural, Urology	16.8 x 23.5 mm FOV 115°	Probe Frequency Range: 1.8 – 8.6 MHz	N/A	Venue, Venue Go and Venue Fit
 12S-RS	Wide-band phased array	Pediatric/Neonatal, Pediatric Cardiac, Neonatal Cephalic	13.2 x 17.6 mm FOV 105°	Probe Frequency range: 2.7 - 12.5 MHz	N/A	Venue, Venue Go and Venue Fit
 8C-RS	Wide-band microconvex array	Abdomen, Adult Cardiac, Adult Cephalic, Gynecology, Musculoskeletal, Neonatal Cephalic, Nerve Block, Pediatric/Neonatal, Pediatric Cardiac, Peripheral Vascular, Thoracic/Pleural, Urology	12 x 22 mm FOV 128°	Probe Frequency Range: 3.1 – 12.9 MHz	N/A	Venue, Venue Go and Venue Fit
 C1-5-RS	Wide-band convex array	Abdomen, Adult Cardiac, Fetal/Obstetrics, Gynecology, Musculoskeletal, Nerve Block, Pediatric/Neonatal, Pediatric Cardiac, Peripheral Vascular, Rheumatology, Thoracic/Pleural, Urology	17.2 x 69.3 mm FOV 59°	Probe Frequency Range: 1.4 – 5.7 MHz	Multi-angle, reusable bracket, disposable sleeve	Venue, Venue Go and Venue Fit
 E8C-RS	Wide-band microconvex array	Gynecology, Obstetrics, Transrectal, Transvaginal, Urology	16.9 x 21.2 mm FOV 128°	Probe Frequency Range: 2.2 – 12.5 MHz	N/A	Venue, Venue Go and Venue Fit
 L4-12t-RS	Wide-band linear array with 4 configurable buttons	Abdomen, Musculoskeletal, Neonatal Cephalic, Nerve Block, Ophthalmic, Pediatric/Neonatal, Peripheral Vascular, Rheumatology, Small Parts, Thoracic/Pleural, Vascular Access	12.7 x 47.1 mm FOV 38.4 mm	Probe Frequency Range: 3.4 – 12.6 MHz	Multi-angle, Free angle and out-of-plane; reusable bracket*, disposable sleeve	Venue, Venue Go and Venue Fit

Transducer	Description	Applications	Footprint	Bandwidth	Biopsy Guide	System Compatibility
 L12n-RS	Wide-band linear array with 4 configurable buttons	Abdomen, Musculoskeletal, Neonatal Cephalic, Nerve Block, Ophthalmic, Pediatric/Neonatal, Peripheral Vascular, Rheumatology, Small Parts, Thoracic/ Pleural, Vascular Access	12.7 x 47.1 mm FOV 38.4 mm	Probe Frequency Range: 3.4 – 12.6 MHz	Multi-angle, Free angle and out-of-plane; reusable bracket*, disposable sleeve	Venue
 12L-RS	Wide-band linear array	Abdomen, Musculoskeletal, Neonatal Cephalic, Nerve Block, Ophthalmic, Pediatric/Neonatal, Peripheral Vascular, Rheumatology, Small Parts, Thoracic/Pleural, Vascular Access	12.7 x 47.1 mm FOV 38.4 mm	Probe Frequency Range: 3.4 – 12.6 MHz	Multi-angle, Free angle and out-of-plane; reusable bracket, disposable sleeve	Venue, Venue Go and Venue Fit
 L4-20t-RS	Wide-band linear array with 4 configurable buttons (XDclear array)	Abdomen, Dermatology, Musculoskeletal, Nerve Block, Neonatal Cephalic, Ophthalmic, Pediatric/Neonatal, Peripheral Vascular, Rheumatology, Small Parts, Thoracic/Pleural, Vascular Access	11.56 x 48.43 mm FOV 38.4 mm	Probe Frequency range: 3.0 - 20 MHz	Multi angle reusable bracket*, disposable sleeve	Venue, Venue Go and Venue Fit
 ML6-15-RS	Wide-band linear array with (active matrix array)	Dermatology, Musculoskeletal, Nerve Block, Neonatal Cephalic, Ophthalmic, Pediatric/Neonatal, Peripheral Vascular, Rheumatology, Small Parts, Thoracic/Pleural, Vascular Access	16.0 x 60.7 mm FOV 50.4 mm	Probe Frequency range: 4.0 - 18.9 MHz	N/A	Venue, Venue Go and Venue Fit
 L8-18i-RS	Wide-band linear array, hockey-stick shape	Abdomen, Dermatology, Intraoperative, Musculoskeletal, Neonatal Cephalic, Nerve Block, Ophthalmic Peripheral Vascular, Pediatric/Neonatal, Rheumatology, Small Parts, Thoracic/Pleural, Vascular Access	11.1 x 34.8 mm FOV 25.2 mm	Probe Frequency Range: 2.5 – 16.8 MHz	N/A	Venue, Venue Go and Venue Fit
 L10-22-RS	Wide-band linear array	Dermatology, Musculoskeletal, Nerve Block, Pediatric/Neonatal, Peripheral Vascular, Rheumatology, Vascular Access	8.0 x 20.3 mm FOV 16.2 mm	Probe Frequency range: 7.0 - 22.2 MHz	N/A	Venue, Venue Go and Venue Fit
 9L-RS	Wide-band linear array	Abdomen, Fetal/Obstetrics, Gynecology, Musculoskeletal, Neonatal Cephalic, Nerve Block, Ophthalmic, Pediatric/Neonatal, Peripheral Vascular, Rheumatology, Small Parts, Thoracic/ Pleural, Urology, Vascular Access	14.1 x 53.0 mm FOV 44.2 mm	Probe Frequency Range: 2.4 – 10.5 MHz	Free angle; reusable bracket, disposable sleeve	Venue, Venue Go and Venue Fit
 6Tc-RS	Wide- band phased array for transesophageal echocardiography	Adult Cardiac	12 x 14 mm with length of 45 mm FOV 90°, Depth of Field 20 cm	Probe Frequency Range: 1.8 – 8.8 MHz	Sheath	Venue, Venue Go and Venue Fit

*Button #1 may not be usable with out of plane needle guide.



© 2021 General Electric Company – All rights reserved.

cSound, Venue, Venue Go and Venue Fit are trademarks of General Electric Company.

Venue Family R3 Release.

Products mentioned in the material may be subject to government regulation and may not be available in all countries. Shipment and the effective sale can only occur after approval from the regulator. Some applications or probes may not be available on all product versions.

Please check your user manual for specific details. Contact your GE representative for details.

DOC2391320